

CS CERAMIC CO.,LTD since 1974

FAX:0086-0731-23685198 TEL:0086-0731-23685198

ADD:Chashan Town Liling City Hunan Province ,China

Consumables list for LECO instrument

FROM:ECHO

Email:info@csceamic.com



Name	Item NO.	Description	PCS/CTN	LECO PART	ELRA PART	ALPHA /HORIBA PART
Carbon sulfur crucible (STANDAR CERAMIC CRUCIBLE)	CS001	C ≦ 3ppm, S ≦ 2ppm, foil-wrapped	1000pcs/CTN	528-018	90150	905.200.380.001
Carbon sulfur crucible (SUPER LOW)	CS002	C ≦ 1ppm, S ≦ 1ppm, foil-wrapped	1000pcs/CTN	528-120		
Crucible Cover	CS003	With 10mm Hole For Use With C4500 & C4501	1000pcs/CTN	528-043 528-040		905.130.200.001
Regular Graphite crucible	CS005	for typical nitrogen/oxygen analyses ,9.86x12.65 mm	1000pcs/CTN	776-247	90190	AR6247/AR053HD
Lower Electrode Tip	CS005-1	for use with 776-247 Crucibles with automation	1/pkg	611-351-181		
Lower Electrode Tip	CS005-2	for use with 776-247 Crucibles without automation	1/pkg	611-351-182		
Graphite crucible	CS005-3	Graphite Crucible 767-277 pack of	1000pcs/CTN	767-277		AR277
Graphite Crucible	CS005-4	Graphite Crucible 760-034 pack of 1000/100	1000pcs/CTN	760-034		
Graphite Crucible (RH-1E)	CSMY-1	RH-1E	1000pcs/CTN	762-213		AR213
Graphite Short Crucible (RH-2)	CSMY-2	RH-2	1000pcs/CTN	769-520		
Graphite Tall Crucible (RH-3)	CSMY-3	RH-2	100pcs/CTN	769-761		

Graphite Cover	CSMY-4	Graphite Cover For Use With C4605 (769-520)	100pcs/CTN	769-569		
Graphite crucible	CSMY-5	760-414 pack of 100/1000	1000pcs/CTN	760-414		AR414
Graphite Crucible	CSMY-6	763-213 pack of 100/1000	1000pcs/CTN	763-213		
Graphite Crucible	CSMY-7	Graphite Crucible 778-674 pack of 100	100pcs/CTN	778-674		
Graphite Crucible RH-1	CSMY-8	769-994 pack of 100	100pcs/CTN	769-994		
Inner Graphite crucible	CS006	OH900,Inner Graphite Crucible for applications requiring uniform crucible with 775-433 Crucibles,10.45*12.65mm	1000pcs/CTN	775-431	90180	AR5431
Outer Graphite Crucible	CS007	outer Graphite Crucible for applications requiring uniform crucible with 775-431 and 775-892 Crucibles,12.85x14.8 mm	1000pcs/CTN	775-433	90185	AR433
Inner Graphite Crucible	CS020	for applications requiring uniform crucible temperature use with 775-433 Crucibles	100/pkg	775-892		
Graphite Crucibles RHEN600-602	CS004	High Temp - for use with EN units (requires 619-896 Tip) 12*15.4	1000pcs/CTN	619-895		
Lower Electrode Tip	CS004-1	for use with the 619-895 Crucible	1/pkg	619-896		
High Temp Graphite Crucible	CS008	OH900/0.5" OD - for high temp alloys and refractories 10.17*12.65mm	1000pcs/CTN	782-720		AR782HD
High Temp Graphite Crucible	CS022	0.625" OD - for high temp alloys and refractories	1000pcs/CTN	782-720S/782-720SHD		AR872HD
Lower Electrode Tip	CS023	for 782-719 and 782-720 Crucibles without automation	1/pkg	782-721		
Graphite Crucible	CSMY23	For LECO OHN-836 Bag of 200	200pcs/CTN	782-719		AR782p

Graphite Crucible	CSMY24	For LECO OHN-836 Bag of 1000	1000pcs/CTN	780-890		AR7808
Graphite Crucible RHN600	CS009	High Temp - for use with EN units (requires 619-896 Tip)	1000pcs/CTN	619-895		AR619
Lower Electrode Tip	CSMY24	for use with the 619-895 Crucible	1/pkg	619-896		
Graphite Crucibles(RH-1EN)	CS010	for use with EN Units (aluminum applications)9.8*12.6mm	1000pcs/CTN	764-330		AR331S
Lower Electrode Tip	CSMY25	for use with 764-330 Crucibles with automation	1/pkg	611-351-180		
Larger ID Crucibles	CS026	typically used for oxygen analysis and in conjunction with 780-892 Covers 1000/pkg	1000/pkg	Larger ID Crucibles		
Graphite Covers/lid	CS027	minimize gettering effect	1000/pkg	780-892		AR0892
High Purity Outer Graphite Crucible	CS028	for oxygen in silicon applications, use with 782-795HD Crucibles	1000/pkg	782-703		AR2703
Inner Crucible	CS029	for oxygen in silicon applications, use with 782-703 Crucibles	1000/pkg	782-795HD		AR2795
Graphite crucible Hydrogen	CSMY-29	Graphite crucible Hydrogen Horiba 905.200.020.001SP pack of 1000	1000pcs/CTN			905.200.020.001SP
Graphite Outer Crucible	CS030	Sold Individually 11.2*14.25mm	1000/pkg			905.200.010.000/AR1533
AR1535 inner Crucible	CS030-1	Suit for 905.200.010.000/AR1533	1000/pkg			
Graphite Crucibles Extra Large	CSMY-30	Horiba EMGA620 905200010001XL	200pcs/CTN			905.200.010.001XL
Graphite crucibles	CSMY-31	H-mat 2020 - S309243000 (Strohlein)	100pcs/CTN			S309243000
Graphite crucibles	CSMY-32	Horiba EMGA620 905200010001	1000pcs/CTN			905.200.010.001
Graphite crucibles	CSMY-33	NH-mat 243 JUWE S309236000 (Strohlein) Pack of 100				S309236000
Inner Graphite Crucible	CSMY-34	Horiba pack of 1000 905.200.050.001				905.200.050.001

Graphite Inner Crucible	CS031	9.88*11.05mm	1000/pkg			905.201.790.001/AR79001
AR31360 Eltra Graphite Tip	CS032	Compatible Instruments OHN2000/ON900	500/CTN		31360	AR31360
Ceramic boat unglazed	CS012	Re-usable ceramic boats/55mm*21mm*11.5mm	100/500pcs/CTN	529-203	90155/90153	
Conbusiton boat	CS013	Low Profile SC144 58mm*23mm*13.5mm	500pcs/CTN	529-204		
Disposable porcelain boats	CS014	86x13x10 mm	500/1000pcs/CTN	528-053	90160	
Horiba Combustion Boat	CS014-1	7.0ml, 97 X 16, Size 6A, 10 mm	1000pcs/CTN			905.201.140.001
Horiba Boat Cover	CS014-2	for No. 6A for C/S Analysis	1000pcs/CTN			905.201.190.001
Fittering crucible	CSMY-15	Filtering Crucible	100PCS/CTN	528-028		
Fittering crucible	CSMY-16	Filtering Crucible	100PCS/CTN	528-030		
Large Crucible	CS018	Large Crucible for TGA 20cc Leco TGA 500 / 501 / 601 / 701, MAC 400 / 500	20PCS/CTN	529-047		AR9047
Large Crucible Cover	CS019	Crucible cover for TGA 20cc	20PCS/CTN	529-048		AR9048
Small Crucible for TGA 16cc	CSMY-17	Small Crucible for TGA 16cc	20PCS/CTN	529-042		AR9042
Small Crucible Cover	CSMY-18	Small Crucible Cover TGA 16cc	20PCS/CTN	529-043		AR9043
Ceramic Crucibles	CS016		50PCS/CTN	529-026		
Porous Crucible	CS017		5PCS/CTN	614-961-110		

REMARKS:

1.Sample:DHL/TNT/UPS/FEDEX Can be chosen ,client need suply freight collect No.

2.Special offer: Most of our crucibles are produced to order, please mark clearly height, bottom diameter, top diameter, wall thickness and any special requirements i.e. lips etc. To speed up quotations please supply a drawing and if appropriate, the name of the system in which they are to be used.(www.cscearmic.com)

CONTACT :info@CSCERAMIC.COM